



# DNABLE<sup>®</sup> Molecular Detection Kit for *Cmm* in Tomato Tissue

Catalog No. DF 091 PT

Part # 12009

## Highlights:

- Molecular detection of *Cmm* in tomato tissue
- Rapid amplification and detection in 10 minute assay

## Contents of DNABLE Kit:

- A. RB1 Reaction Buffer
- B. *Cmm* Master Mix
- C. Flat Caps
- D. MB1 Extraction Buffer
- E. 1.5 mL Extraction Tubes
- F. Blue pestles

## Materials Not Provided:

- 3 mm Harris punch\* (Part# 11209)
- Minipets: 5  $\mu$ L and 50  $\mu$ L\* (Part# 11818, 11203)
- Marker
- 8-well AmpliFire Reader\* (Part# 11974)

\*available through EnviroLogix

## Intended Use

This test kit is intended for qualitative detection of *Clavibacter michiganensis michiganensis* (*Cmm*). The results of this test may facilitate rapid, point of need detection of *Cmm* in tomato plants.

## How the Test Works

DNABLE is an isothermal nucleic acid amplification technology enabling rapid amplification of a specific DNA target. In this test, samples are collected, processed, and added to the reaction buffer. The reaction buffer containing sample is then transferred to the lyophilized master mix, containing all the reagents needed to specifically recognize, amplify and detect the *Cmm*-specific DNA.

The amplified *Cmm*-specific DNA is detected in real-time and the results are displayed and interpreted within 10 minutes using our 8-well DNABLE Reader.

## Precautions and Notes

DNABLE is a highly sensitive assay. Therefore the following precautions are recommended to reduce the chance of sample contamination:

- Clean the work stations and pipettes before and after use with 10% bleach solution
- It is recommended to physically separate sample preparation activities from DNABLE assay activity
- Do not reuse kit disposables
- Use fresh pipette tips for each sample
- Discard used tips in a sealed container containing 10% bleach solution
- Use careful pipetting techniques to avoid cross-contamination between samples; avoid reaching over or pipetting over open tubes
- Wear disposable gloves and change between handling of samples

**Important: Never open reaction tubes after reaction has occurred, as this will release amplified material into the environment and may contaminate**

subsequent reactions. Care should be taken when disposing of run reaction tubes to avoid possibility of tube leakage. Place completed reaction tubes back in original zippered pouch prior to disposal.

## Kit Components

- A. RB1 Reaction Buffer: Provided in green 8-well strip tubes (3 strips)
- B. *Cmm* Master Mix: Lyophilized reagents provided in clear 8-well strip tubes (3 strips)
- C. Flat Caps: used for capping the clear tubes prior to assay start (3 strips)
- D. MB1 Extraction Buffer: One dropper bottle containing 15 mL of extraction buffer (for sample preparation)
- E. 1.5 mL Extraction Tubes (25): One bag of 25 tubes for sample extraction
- F. Blue Pestles (25): One bag of 25 pestles for sample extraction



## Before Testing

- Remove needed DNAble Kit reagents from refrigerated storage. Allow reagents to come to room temperature before opening sealed white pouches.
- Ensure that all assay reagents, extracted sample, pipettes and flat caps are ready for use.
- Wear gloves and wipe gloves with disinfecting wipes between tomato tissue sample collection.



## Tomato Tissue Sample Collection and Preparation

1. Sandwich a section of stem or leaf tissue between the cap of a closed extraction tube and the Harris Uni-Core punch and apply pressure to the stem while rotating the Uni-Core back and forth.

Alternatively, harvest a section of stem, slice open lengthwise, and remove sample from tissue inside stem using Uni-Core punch as above.

*Note: Like all pathogen testing, sampling symptomatic tissue increases the likelihood of detecting Cmm with DNAble.*

2. Release the cored sample in the 1.5 mL extraction tube by depressing the spring action plunger. Verify that the sample is at the bottom.

*Note: Clean the cutting edge of the Uni-Core punch between each sample: dip into 10% bleach solution (or other suitable sanitization solution), depress the plunger several times, then rinse with potable water.*

3. Insert the pestle into the tube and grind the tissue by rotating the pestle against the sides of the tube with twisting motions. Continue this process for 5 seconds or until the plant tissue is gently ground to release bacteria.
4. Holding the dropper bottle vertically, add 0.5 mL (to the 0.5 line on tube) of MB1 Extraction Buffer into the tube containing the tissue sample. Do not touch the dropper or pipette tip to the tube or tissue.
5. Repeat grinding step briefly to mix tissue with extraction buffer. Dispose of pestle.
6. Using **red MiniPet** (or pipette set to 5  $\mu$ L), **add 5  $\mu$ L of sample extract to 1 well of the green reaction buffer strip**. Repeat with remaining samples and wells, changing pipette tips between samples.




Take tissue



Release tissue in tube



Grind tissue gently



Add buffer, grind again briefly to mix.

**Pipette 5  $\mu$ L sample in one well**  
**Change pipette tips between samples/wells**





## How to Run the DNABLE Assay

- On the AmpliFire screen, select “Execute Reaction” then select “Scan Product Code”. Use the barcode on the master mix foil pouch to scan the Cmm protocol on the 8-well AmpliFire Reader. Cmm\_Lot # will display. Select “Next”.
- Under “Reaction Name” enter an appropriate reaction description. **This description is placed at the beginning of the file name.** Select “Next”.
- To enter sample specific information, add sample descriptions to the screens for Wells 1 through 8, clicking “Next” to advance to each Well. Select “Finish” to skip well-specific sample entry.
- Remove clear Master Mix tubes from the foil pouch and gently tap down to ensure that the white pellet is at the bottom of the tubes.
 

*Important: Mark flat cap for orientation of the clear Master Mix tubes (writing on clear tubes will interfere with results interpretation or leave marker residue in instrument).*
- Using the **yellow MiniPet** or a multichannel pipette set to **50  $\mu$ L**, transfer **50  $\mu$ L from green strip tubes (containing sample)** to clear Master Mix tubes. Do **not** mix within the clear tube.
- Cap** Master Mix tubes with provided **Flat Caps** strip.
 

*Important: Ensure that the tubes are completely sealed with flat caps*
- Gently flick down on the resuspended, capped master mix to ensure that no bubbles are at the bottom of the tube and that master mix is fully resuspended.
- Inspect tube to ensure that **no air bubbles are present within the sample volume** (a bubble at the top is fine) and that **cap is completely sealed**.
- When the strip is ready select “Start”. Place resuspended, capped clear strip tube in instrument and press “Ok”.
- After 10 minutes, the AmpliFire will produce a short beeping sound and display final results. Results will be interpreted as Not Detected (-) or Positive (+).
 

*Important: Positive results may be interpreted prior to assay completion, but the full assay time must be complete for complete result interpretation. (Empty wells will be interpreted as negative.)*
- After the assay is complete, carefully **remove run reaction strip tubes from instrument and place in opened foil pouch** (used to store master mix), seal and discard in waste container.
- To export results, return to the home screen, then “View Results”. Insert a USB storage device into instrument (left side) and select each run to export and “Export Selected” and “OK.” The results will be saved in a PDF summary report as well as .csv file format.



**LIMITED WARRANTY**

EnviroLogix Inc. (“EnviroLogix”) warrants the products sold hereunder (“the Products”) against defects in materials and workmanship when used in accordance with the applicable instructions for a period not to extend beyond a product’s printed expiration date. If the Products do not conform to this Limited Warranty and the customer notifies EnviroLogix in writing of such defects during the warranty period, including an offer by the customer to return the Products to EnviroLogix for evaluation, EnviroLogix will repair or replace, at its option, any product or part thereof that proves defective in materials or workmanship within the warranty period.

ENVIROLOGIX MAKES NO OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. The warranty provided herein and the data, specifications and descriptions of EnviroLogix products appearing in EnviroLogix published catalogues and product literatures are EnviroLogix’ sole representations concerning the Products and warranty. No other statements or representations, written or oral, by EnviroLogix’ employees, agents or representatives, except written statements signed by a duly authorized officer of EnviroLogix Inc., are authorized; they should not be relied upon by the customer and are not a part of the contract of sale or of this warranty.

EnviroLogix does not warrant against damages or defects arising in shipping or handling, or out of accident or improper or abnormal use of the Products; against defects in products or components not manufactured by EnviroLogix, or against damages resulting from such non-EnviroLogix made products or components. EnviroLogix passes on to customer the warranty it received (if any) from the maker thereof of such non-EnviroLogix made products or components. This warranty also does not apply to Products to which changes or modifications have been made or attempted by persons other than pursuant to written authorization by EnviroLogix.

**THIS WARRANTY IS EXCLUSIVE.** The sole and exclusive obligation of EnviroLogix shall be to repair or replace the defective Products in the manner and for the period provided above. EnviroLogix shall not have any other obligation with respect to the Products or any part thereof, whether based on contract, tort, strict liability or otherwise. Under no circumstances, whether based on this Limited Warranty or otherwise, shall EnviroLogix be liable for incidental, special, or consequential damages.

This Limited Warranty states the entire obligation of EnviroLogix with respect to the Products. If any part of this Limited Warranty is determined to be void or illegal, the remainder shall remain in full force and effect.



For Technical Support  
Contact Us At:

**EnviroLogix Inc.**  
500 Riverside Industrial Parkway  
Portland, ME 04103-1486 USA  
Tel: (207) 797-0300  
Toll Free: 866-408-4597  
Fax: (207) 797-7533

e-mail: [dnable@envirologix.com](mailto:dnable@envirologix.com)

Website: [www.envirologix.com](http://www.envirologix.com)

**LICENSE**

EnviroLogix has developed this kit using proprietary technologies licensed from Ionian Technologies, and the University of Medicine and Dentistry of New Jersey. Dye compounds in this product are sold under license from Biosearch Technologies, Inc. and protected by U.S. and world-wide patents either issued or in application. The license covers agricultural (including GMO), plant pathogen and veterinary use only and may not be used for human *in vitro* diagnostic (IVD) applications.

*EnviroLogix, the EnviroLogix logo, DNable, and the DNable logo are trademarks of EnviroLogix Inc.*

© EnviroLogix 2015